

# STEREO DISSECTING MICROSCOPE

## ABRON MODEL AB-31A

### Sliding Objectives



- \* Three Dimensional Stereo & Sharp Viewing.
- \* Useful for Cut Flowers, Plant Parts / Whole Specimens.
- \* Binocular Head with Inter Pupillary Distance 55-75mm with Eye Diopters Adjustment.
- \* different Magnification , Superior Optics and Prisms.

Eye piece	10x (paried)		
OBJECTIVE	2X	Magnification	20x
OBJECTIVE	4X	Magnification	40x

# STEREO DISSECTING MICROSCOPES

## MODEL AB-31B

### Sub Satage Light variable

- \* Three Dimensional Stereo & Sharp Viewing.
- \* Useful for Cut Flowers, Plant Parts / Whole Specimens.
- \* Binocular Head Prismaic with Inter Pupillary Distance 55-75mm with Eye Diopters Adjustment.
- \* different Magnification , Superior Optics and Prisms.

Eye piece	10x (paried)		
OBJECTIVE	2X	Magnification	20x
OBJECTIVE	4X	Magnification	40x

