



**AB-34T STEREO ZOOM TRINOCULAR MICROSCOPE**

Trinocular Head	Angle of viewing:45°
Zoom ratio:	1:6.4
Eye Piece	WF10X/20mm
Objective:	0.7x-4.5x
Magnification	7x-45x
Interpupillary Distance	adjustment range:55-75mm
Working distance:	100mm
Reflected transmitted	LED illumination

- |           |                    |        |              |            |               |
|-----------|--------------------|--------|--------------|------------|---------------|
| AB-34B-D1 | Stereo with Camera | AB-14a | <b>1.5MP</b> | + software | CD for laptop |
| AB-34B-D3 | Stereo with Camera | AB-14a | <b>3MP</b>   | + software | CD for laptop |
| AB-34B-D5 | Stereo with Camera | AB-14a | <b>5MP</b>   | + software | CD for laptop |
| AB-34B-D6 | Stereo with Camera | AB-14a | <b>5MP</b>   | + software | CD for laptop |
- HIGH RESOLUTION CAMERA**