

A b

**Conductivity pH mV Digital LED Temperature Meter Table with combined and other electrode**

**pH Range x Accuracy** 0 to 14 pH x 0.01 pH  
**Temp Compensation** Automatic with Manual  
**Buffer recognition** 4.00, 7.00, 9.20 auto  
**Slope adjustment** 80% to 120%  
**Calibration Auto/ manual** 2-point,3-point

**mV (ORP) Range** 0 to ± 1999 mV X ±1mV  
**Calibration** manual

**Temperature Range** -20 to +200°C x 0.1' C  
**Input type Pt 100** Calibration manual  
**Power Supply** 230 V ± 10% AC 50Hz

**Conductivity Mhos:** Mhos/Cm To 20m  
**Accuracy** 0.3% F.S. 2 Micro Mhos/Cm To 20m.Mhos/Cm 1% Fs In Last Range 20m Mhos/Cm  
**Resolution** 0.1 Micro Mhos/Cm  
**Measuring Frequency** 1000hz  
**Temperature Compensation** 0°C To 50°C Manual  
**Cell Constant** 0.4 To 1.6 Adjustable  
**Display** 3½ Digit Led Display

- Accessories**
- 1) Ph Electrode
  - 2) Conductivity Cell
  - 3) Temperature Probe
  - 4) Electrode Stand
  - 5) Dust Cover
  - 6) Instruction Manual

Cat No	Specifications	Range	Accuracy
AC-308-DC	pH meter mV temp Conductivity table model with combined electrode	0-14pH	0.01pH
AC-308-DP	pH meter mV temp Conductivity+ printer output table model with combined electrode	0-14pH	0.01pH
AC-308-E	pH meter Spare combined electrode	0-14pH	0.01pH



**Conductivity pH mV ORP TDS Temperature Meter Touch Table alpha numeric display combined electrode Printer Out RS232**

**Display Format** simultaneously auto Setting  
**Storage** 1000 samples  
**Interfacer** RS-232, Real Time Display  
**Error Indication** pH  
**pH Range x Accuracy** 0 to 14 pH x 0.01 pH  
**Temp Compensation** Automatic with Manual  
**Buffer recognition** 4.00, 7.00, 9.20 auto  
**Slope adjustment** 80% to 120%  
**Calibration Auto/ manual** 2-point,3-point  
**mV (ORP) Range** 0 to ± 1999 mV X ±1mV  
**Calibration** manual  
**Temperature Range** -20 to +200°C x 0.1' C  
**Input type Pt 100** Calibration manual  
**Power Supply** 230 V ± 10% AC 50Hz  
**Conductivity Mhos:** Mhos/Cm To 20m x 2 Mhow  
**Last Range** 20m Mhos/Cm x 0.1  
**Measuring Frequency** 1000hz  
**Temperature Compensation** 0°C To 50°C Manual  
**Cell Constant** 0.4 To 1.6 Adjustable

**Accessories**

- 1) pH Electrode-mV
- 2) Standard Solution
- 3) Temperature Probe
- 4) Electrode Stand
- 5) Dust Cover
- Conductivity Cell TDS 7.0 pH & 4.0 pH
- 6) Instruction Manual

Cat No	Specifications	Range	Accuracy
AC-308-NP	pH meter mV ORP temp micro + printer output table model with combined electrode	0-14pH	0.01pH
AC-308-N	pH meter mV ORP temp micro no printer output table model with combined electrode	0-14pH	0.01pH



**pH Meter Abron Spares**

Cat No	Specifications	Range	Accuracy
AC-307-SP	pH meter spare standard solution sachet	4,7,9.2pH	
AC-307-SB	pH meter bottles with capsule	4pH 7pH	
AC-308-E	pH meter Spare combined electrode	0-14pH	0.01pH



ABRON EXPORTS  
www.abronexports.com

**pH Meter Abron Spares**

Cat No	Specifications	Range	Accuracy
AC-307-SP	pH meter spare standard solution sachet	4,7,9.2pH	
AC-307-SB	pH meter bottles with capsule	4pH 7pH	
AC-308-E	pH meter Spare combined electrode	0-14pH	0.01pH



ABRON EXPORTS  
www.abronexports.com