





Abbron's

ISO 9001-2008
An Certified Company

NEW PRODUCTS

Abbron
Instruments

BINOCLUAR MICROSCOPE



Model AB-12 A

Specification: Binocular Research Microscopes

- | | |
|--|---|
| • 45mm DIN Achromatic Objective | : 4X,10X,40X(S),100X(S, Oil) |
| • Wide Field Plane-Scope Eyepiece | : WF10X |
| • Total Magnification | : 40X-1000X |
| • Articulated Free Binocular Head, Interpupillary Distance | : 55-75mm, 30° Inclined, 360° Rotatable |
| • Double Layers Mechanical Stage | : 125X125mm, Moving Range X-Y: 60X40mm |
| • Ø2-Ø30mm Iris Diaphragm And Ø32 Filter, ABBE Condenser | : N.A. 1.25 |
| • Coaxial Coarse And Fine Focusing Adjustable Mechanism | : 14mm |
| • Built-in Adjustable Brightness Halogen Lamp
(or Incandescent or Fluorescent or LED) | : 6V/20W |

Price Rs. 25000/-

Abbron Exports, INDIA