

**Oil Bath Digital PID Abron High Temperature Water Bath 350°C x1% Electric with stirrer**

Outer Body Double Walled Powder Coated  
Inner Chamber Made Of Stainless Steel Sheet  
**Heating** By Immersion Type Heater  
220v 50hz Ac 2 Meter 3 Core with Plug  
**Digital Temperature** 5° C Above Ambient 350° C  
**Accuracy** 2-5°c Digital PID Least Count 1°c  
**Stirrer** with Variable Speed

Cat No	Specifications	Type	Size
AC-379-D10	Oil bath Electric Digital PIC Chamber	SS	10"X10"
AC-379-D17	Oil bath Electric Digital PIC Chamber	SS	17"X10"

**Optional Spares**

AC-379-S1	Oil Bath spare silicon oil 1 litre	2x500 ml
AC-379-S5	Oil Bath spare silicon oil 5 litre	5000 ml



**Abron Exports**  
www.abronexports.com

**Oil Bath Digital PID Abron 200°C Temperature Bath x1% with window Electric with stirrer**

Outer Body Double Walled Powder Coated  
Inner Chamber Made Of Stainless Steel Sheet  
**Window** Glass to view chamber  
**Heating** By Immersion Type Heater  
220v 50hz Ac 2 Meter 3 Core with Plug  
**Digital Temperature** 5° C Above Ambient 200° C  
**Accuracy** 2-5°c Digital PID Least Count 1°c  
**Stirrer** with Variable Speed

Cat No	Specifications	Chamber	Size
AC-379-D10W	Oil bath Electric Digital PIC Window	SS	10"X10"
AC-379-D17W	Oil bath Electric Digital PIC Window	SS	17"X10"



**Abron Exports**  
www.abronexports.com

**Oil Bath Thermostatic Control Round High Temperature Bath 350°C Insulated Electric**

Outer Body Double Walled Powder Coated  
Inner Chamber Made Of Stainless Steel Sheet  
**Heating** By Immersion Type Heater  
220v 50hz Ac 2 Meter 3 Core with Plug  
**Temperature-Thermostatic** 5° C Above Ambient 350° C  
**Accuracy** 5°c

Cat No	Specifications	Chamber	Dia x H
AC-379-R10	Oil bath Electric Thermostatic 350°C	SS	10"X10"
AC-379-R12	Oil bath Electric Thermostatic 350°C	SS	12"X12"
AC-379-R14	Oil bath Electric Thermostatic 350°C	SS	14"X12"
AC-379-R18	Oil bath Electric Thermostatic 350°C	SS	18"X12"

**Optional Spares**

AC-379-S1	Oil Bath spare silicon oil 1 litre	2x500 ml
AC-379-S5	Oil Bath spare silicon oil 5 litre	5000 ml



**Oil Bath Abron Digital PID Controlled Round High Temperature Bath 350°C Electric**

Outer Body Double Walled Powder Coated  
Inner Chamber Made Of Stainless Steel Sheet  
**Heating** By Immersion Type Heater  
220v 50hz Ac 2 Meter 3 Core with Plug  
**Temperature-Digital PID** 5° C Above Ambient 350° C  
**Least Count** 1°C

Cat No	Specifications	Chamber	Dia xH
AC-379-R10D	Oil bath Electric Digital PID 350°C	SS	10"X10"
AC-379-R12D	Oil bath Electric Digital PID 350°C	SS	12"X10"
AC-379-R14D	Oil bath Electric Digital PID 350°C	SS	14"X10"
AC-379-R18D	Oil bath Electric Digital PID 350°C	SS	18"X10"

**Optional Spares**

AC-379-S1	Oil Bath spare silicon oil 1 litre	2x500 ml
AC-379-S5	Oil Bath spare silicon oil 5 litre	5000 ml



**Abron Exports**  
www.abronexports.com

**Abron Exports**  
www.abronexports.com