

# Micropipette Digital Variable Volume Special Two Button Aircushion Micropipette Abron AC-289SV



- Variable with Eject Tips Easy setting by simply turning mechanical
- Click Stop System for Aqueous Liquids Two Button operation
- Easy to calibrate and service. Pipetting of aqueous liquids,
- two button operation Volume range as below
- Available options: 1-channel Autoclavable for lab uses &
- Purity grades of consumables: PCR clean and DNA forensic
- Autoclavable tip cone and lower section
- Accurate \* Conform to DIN 12650 Standards

Cat No	Type	Volume ml	Least count LC
AC-289-SV10	Micro Pipette Digital variable	0.5-10 ul	0.1 ul
AC-289-SV20	Micro Pipette Digital variable	2-20 ul	0.1 ul
AC-289-SV50	Micro Pipette Digital variable	5-50 ul	1 ul
AC-289-SV1H	Micro Pipette Digital variable	20-100 ul	1 ul
AC-289-SV2H	Micro Pipette Digital variable	20-200 ul	1 ul
AC-289-SV1T	Micro Pipette Digital variable	200-1000 ul	10 ul
AC-289-SV1TA	Micro Pipette Digital variable	100-1000 ul	10 ul
AC-289-SV5T	Micro Pipette Digital variable	1-5 ml	100 ul
AC-289-SV10T	Micro Pipette Digital variable	2-10ml	100 ul