

Table top Autoclave front loading 23Ltr Class B Front Loading Dental Autoclave Abron AC-325TB



Features

- B type front loading autoclave with triple Vacuum cycle.
- Friendly operation interface, simple and easy.
- B & D test to check temperature penetration capacity of sterilizer Leak test to monitor whether there is leakage.
- 3 temperature selection : 121° C , 134° C & B+
- Sterilization of unwrapped/wrapped, solid/hollow, porous/non-porous instruments can be done.
- Step button to select a date, month, year and time to store the sterilization cycle.
- 2 separate water tanks available: Distilled water with 3.5litre capacity and waste-water tank with 5litre capacity.
- Cooling condenser in the waste-water tank to reduce the water temperature and prevent deformation of the tank.
- Water quality check sensor in the distilled water tank available.
- An optional Printer is available.
- Double door lock system i.e. airlock & mechanical lock to prevent potential damage.
- B+ cycle in 134° C temperature and sterilization holding time of 18 min for effective sterilization of Kits.