

Microscope Olympus

Microscope Student ISI High Quality Optics

AV-16-ISI

Huygenian Eye Piece	Achromatic Objective	Magnification
10X	10 X	100 X, 150x
15X	45 X	450X , 675X

Viewing Head : Monocular 90°
Nosepiece : Triple revolving Upright
Objective : Achromatic 10x, 40x Spring Loaded SL
Stage : 120x120mm with stage clip
Focusing Knobs : Coarse and Fine separate
Condenser : Iris Diaphragm
Illumination : Plano-concave mirror
Carrying Case Styrofoam **Supplied with:** Dust cover

Cat No **Type** **Optics**
 AB-16-ISI Microscope student Iris diaphragm high quality lens
 Eye piece 10x,15x Objective 10x,45x



Abron Exports
www.abronexports.com



Abron Exports
www.abronexports.com

Optional Accessories:

AB-57-C LED Lamp, Variable 2x1.5v Cell
 AB-54-C Eye Piece 10X WF

Student Microscope Olympus HSA

Huygenian Eye Piece	Achromatic Objective	Magnification
10X	10 X	100 X, 150x
15X	45 X	450X , 675X

Viewing Head : Monocular 90°
Nosepiece : Triple revolving Upright
Objective : Achromatic 10x, 40x Spring Loaded SL
Stage : 120x120mm with stage clip
Focusing Knobs : Coarse and Fine separate
Condenser : Iris Diaphragm
Illumination : Plano-concave mirror
Carrying Case Styrofoam **Supplied with:** Dust cover

Cat No **Type** **Optics**
 AB-03-A Microscope student Iris diaphragm high quality lens
 Eye piece 10x,15x Objective 10x,45x



HSA

Student Microscope bron

same as

Olympus HAS AB-03-B

Huygenian Eye Piece	Achromatic Objective	Magnification
10X	10 X	100 X, 150x
15X	45 X SL	450X , 675X

Viewing Head : Monocular tube
Nosepiece : Triple revolving Upright
Mechanical Stage : 120x120mm with stage clip
Focusing Knobs : Coarse and Fine separate
Condenser : Fixed Abbe N A 1.25 with Iris Diaphragm
Illumination : Plano-concave mirror
Carrying Case : Styrofoam Thermocole
Supplied with : Manual & Dust cover

Cat No **Type** **Optics**
 AB-03-B Microscope student Abron same as HAS
 Iris diaphragm high quality lens
 Eye piece 10x,15x Objective 10x,45x



Abron Exports
www.abronexports.com

Single Nose Abron Microscope AB-09-A

Cat No **Type** **Optics**
 AB-09-A Single nose Microscope Eye piece 10 X objective 10 X

Huygenian Eye Piece	Achromatic Objective	Magnification
10X	10 X	100 X

Nosepiece : Single tube
Mechanical Stage : 90 x 90mm with stage clip
Focusing Knob : Coarse and Fine separate
Condenser : Fixed Abbe N A 1.25 with Iris Diaphragm
Illumination : Plano-concave mirror

