

Inclined Illuminator Microscope

**360° Rotatable Head
Student
Model AB-02B**

STUDENT MICROSCOPE

**360° ROTATABLE HEAD
ABBA CONDENSER LENS
WITH IRIS DIAPHRAGM**

- Triple Revalving nose
Fixed condenser
With Iris Diaphram
- Sturdy body,
Horse Shoe base.
- Coarse
& Fine focussing knobs
- Stage clips,
Reflector,
Dust cover, Cloth.

HUYGENIEN EYE PIECE
10x and 15x

ACHROMATIC OBJECTIVES
10x and 45x

MAGNIFICATIONS
100X 150X
450X, 675X

