

Abron Kymograph Part Smoke Burner & Holder Stand Setup AB-2550BRD



This is useful instrument Set for Paper Blackening
For Arrow Making like graph on paper on Pharmacy Study

Parts

Smoke Burner With Handle to lift Wick Up Down
Die cast powder coated
with NP rod for drum
and stand

For Kymograph experiment

<u>Cat No</u>	<u>Specification</u>	<u>Type</u>	<u>Unit</u>
AB-2550BRD	Kymograph smoke burner with Drum holder		
AB-2550BR	Kymograph smoking burner for Drum		
AB-2550BRW	Kymograph smoking burner Wick SPARE		