



**BE065 Basic Microscope, Abron**  
**40 x 100 x 400**

all metal body,  
adjustable stop to prevent damage slide & lens  
triple revolving nose piece with piece  
with positive stop and  
4x, 10x and 40x Acromatic objectives.  
10 x wide field eye piece  
88 x 82 mm metal stage  
with lock on stage clips  
& rotary disc diaphragm,  
42mm planoconcave mirror illuminator  
with inclined eye tube

- Sturdy body, Horse Shoe base.
- Coarse & Fine focussing knobs.
- Stage clips, Reflector, Dust cover, Cloth.

**Huygenion Eye Piece** 10x WF

**Acromatic Objectives** 4x 10x 40x

**Magnification Range** 40x 100x 400x

**Condenser** Fitted with NA 1.2

**Rotary Disc Diaphragm**

**ABRON EXPORTS**  
A SCIENTIFIC INSTRUMENTS FACTORY



info@abronexports.com  
www.abronexports.com  
AMBALA CANTT-133001