

Haemocytometer Blood counting chamber Set

Neubauer set Bright lined, counting chamber, cover glasses, one red cell pipette, one white blood cell pipette, tubing, mouthpieces and carrying case. Counting leucocytes, erythrocytes and thrombocytes) also count bacteria and fungus spores.

Cat No	Item	Type	Unit
AB-75-HCN	Haemocytometer Neubauer chamber complete	1 set	
AB-75-HCG	Haemocytometer Neubauer chamber German	1 set	
AB-75-HCS	Haemocytometer Cover slip 22x40mm	37Pc	1 box
AB-75-HWB	Haemocytometer Spare RBC Pipette with Tube	1 Pc	
AB-75-HWB	Haemocytometer Spare WBC Pipette with Tube	1 pc	



Haemoglobinometer Sahli's Type complete set

The two side windows are fitted with colour standard strips and middle window is meant for the square tube for colour comparison. The size of all the three windows are Height 3cm and Width 0.5 cm. (Minimum). Square H.B. Tube: 2 No

Cat No	Item	Type	Set
AB-75-HG	Haemoglobinometer with Pipette complete		set
AB-75-HG	Haemoglobinometer Indian Sahli's Type		set
AB-75-HP	Haemoglobinometer Spare Tube WBC+ RBC		set



Hemoglobinometer Hb Strip type Digital AC-309P

Accurate 10 µL capillary venous blood required.
Wide Hb measurement Range 4.5-25.6 g/dL
Hct range of 13-75% Stores up to 1,000 results
battery or optional AC adapter
Wavelengths 525 nm PC Interface
Calibration Automatic Power Source 3 AAA Batteries

Cat No	Item	Type	Set
AC-309P	Haemoglobinometer Digital 10 Strip complete		set



Hemoglobinometer Digital LCD AC-309H 16x2 lines

Microprocessor Hemoglobin Meter Battery Operated with 16x2 lines LCD display for direct HB results soft touch keys, auto calibration. Standard memory Facility 1 ml Sample measurement

Measuring method	Cyanmethemoglobin
Measuring range	0-40 gm/dl of Hb
Dilution ratio	1.25 Sample Vol 1 ml
Display	16x2 Lines LCD Display
Wavelength	530 nm
Keyboard	Soft Touch 4 Keys
Detector	Sensitive Silicon Photodiode
Power Source	230V ± 10% 50 Hz AC
Accessories	1) Set of 5 Test Tube 2) Tube Stand 3) Spare Lamp 4) Dust Cover 5) Instruction Manual.

Cat No	Item	Type	Set
AC-309H	Hemoglobinometer digital Microprocessor LDC 2 line x 16 ch display 0-40 complete		1 set



Hemoglobinometer Microprocessor

LCD display for direct HB results soft touch keys auto calibration. Standard memory Facility 1ml Sample measurement

Measuring method	Cyanmethemoglobin
Measuring range	0-40 gm/dl of Hb
Dilution ratio	1.25 Sample Vol 1 ml
Display	16x2 Lines LCD Display
Wavelength	530 nm
Keyboard	Soft Touch 4 Keys
Zero Setting	Automatic Calibration Auto
Detector	Sensitive Silicon Photodiode
Power Source	230V ± 10% 50 Hz AC
Accessories	1) Set of 5 Test Tube 2) Tube Stand 3) Spare Lamp 4) Dust Cover 5) Instruction Manual.

Cat No	Item	Type	Set
AC-309HM	Hemoglobinometer digital Microprocessor 2 line x 16 ch display 0-40 complete		1 set



Abron Exports
www.abronexports.com