

GM Counter (Geiger- Muller) 6 Digital setup Timer, GM End Tube, stand, Beta Source, stand Absorbers AE-1280A

For counting radiations (e.g. from radioactive substances)
 EHT Power Supply 150 Volts to 1500 Volts,
 Counting System with timer Digital LED
 GM Tube Stand for GM Tube
 Beta Source Absorber plates set
Time Indicator indicates period which scaler accepts
Switch Reset counter system to set zero Fresh reading
Power Need 230 Volts AC 50 Hz
 To detect Alpha particulars, Beta particulars, Gamma
Source supplier Beta as available for research use
 Radiation Source supplied by Bhaba Atomic Research
 Centre against college requirements directly. However, at time
 it can be supplied if available in small Qty special container lid
Warning Keep Radiation Source away from children
Complete setup with Manual Easy to setup & USA

AE-1280P To study The absorption of beta rays and hence
 Determine half-life Value thickness & absorption
 Coefficient of plastic GM Counter setup
 AE-1280R To study statics in Radioactive Measurement
 Half life time GM Counter setup



Cat No	Item	with Power Supply	meters
AE-1280-A	GM Counter (Girger Counting System) 6 digital		
	G.M. Counter with Built in Power Supply setup		
AE-1280-B	GM End Window G.M. Detector		
AE-1280-C	GM Beta sources		
AE-1280-D	GM Gamma reference sources Not now		
AE-1280-E	Aluminium / Lead absorber disc set		



Cat No	Item	with Power Supply	meters
AE-1280-F	Lead sheet for Body		
AE-1280-H	GM Counter Hand Held Portable		

