

# Projection Microscope 15cm screen with Light AB-14P

## AB-14p PROJECTION MICROSCOPE STUDENTS

**Abron**

Ideal for group  
teaching,

**Focussing-** for coarse/fine

**Mechanical Stage-**145x 125mm

**Monocular Head-High Resoulution**

**Lighting :** A built-in base

Halogen 6V20w supply220v AC

**EYE PIECES-** Wide field 10x

**OBJECTIVE-** 10x 20x,

**OPTIONAL OBJECTIVE-**

10x 20x,

**Magnification** 100x

200x

MODELS

AB-14p SIMPLE STAGE

AB-14R MECHENICAL STAGE

