

A b

**Colony Counter Digital 1999 Table Model**

for rapid and accurate counting of bacterial colonies with 4 Digit Display (999max) Dish Size 110mm Dia and Magnification and Pen for counting

**AC-302D4 Colony Counter table Model Table Model with special pen and magnifier with stand**


Digital Display	4 digit (999max)
Dish Size	110mm Dia
Magnification	X1.7
Power	230V±10PHz AC
Size	260x255x170mm
Accessories	1) Marking pen 2) Magnifier Lens with Stand 3) Lens holder 4) Straight Rod 5) Dust Cover 6) Instruction Manual

**Colony Counter Digital 999999 Table Model microprocessor based**

for rapid and accurate counting of bacterial colonies with 4 Digit Display (999max) Dish Size 110mm Dia and Magnification and Pen for counting

**AC-302D6 Colony Counter digital 999999 Table Model with special pen and magnifier with stand**

Digital Display	6 digit (999999max)
Dish Size	110mm Dia
Magnification	X1.7
Power	230V±10PHz AC
Size	260x255x170mm
Accessories	1) Marking pen 2) Magnifier Lens with Stand 3) Lens holder 4) Straight Rod 5) Dust Cover 6) Instruction Manual



**Spares**

- AC-302-P Colony Counter spare Pan
- AC-302-M Colony Counter magnifier