

MICROTOMES Histopathology

Microtome Hand with 10 micron steps + Razer

Sample holder clamp with 10 micron moving by screw
With blade holder razer Table clamp model,
cutting sections of unembedded as well as embedded.
with table clamp device razor Plano concave.


Cat No	Item	Type	Thin	Rs
AB-87-A	Microtome Hand Operated Steps		10 micron	
AB-87-AT	Microtome Hand Operated Steps With table clamp		10 micron	



Microtome Abron Student Rocking Open

heavy base Manual feeding by screw Cambridge type
Suitable for comparatively soft tissue embedded
in paraffin wax sections are uniform in thickness
Salient features: Feed adjusted up to 2 to 20 micron
In step of 2 micron total excursion 20 mm
one razor, three block holder one sharpening

Cat No	Item	Type	Thin
AB-88-A	Microtome Student Open rocking Manual 2 um Feed With Razer , 3 block holder sharpening stone		
AB-88-A	Microtome Student Open rocking Manual 2 um Feed With Blade Razer , 3 block holder sharpening stone		



Microtome Knife with handle Abron

Spare knife suitable for section cutting
Now disposable knives for soft tissue / hard tissue also with us

Cat No	Item	Type	Size	Rs
AB-87-B	Microtome Knife with handle		120mm	
AB-87-C	Microtome Knife with handle		140mm	
AB-87-D	Microtome Knife with handle		160mm	
AB-87-E	Microtome Knife with handle		180mm	
AB-87-F	Microtome Knife with handle		200mm	



Microtome Large Rotary feed 25x2 micron

Closed ERMA Type for cutting of uniform sections up
to 25 microns in steps of 1 micron each The internal structure
is completely encased in a metallic cover selected alloys are
used in fabrication of the machine
With Razor 120mm long, Sharpening stone Oil
Object holder Wheel to cut sections

Cat No	Item	Type	Thin	Rs
AB-89-A	Microtome Rotary close pattern Feeding 2-25 micron 2 micron steps			



Microtome Senior Rocking 0-25 x 1micron

Cat No	Item	Type	Thin	Rs
AB-89-B	Microtome Senior feeding 1-25 micron 1 Micron With Razor 120mm long, Sharpening stone Oil Bottle 60ml, Object holder Wheel to cut sections			



Microtome Automatic Razor Sharpener

for sharpening of all microtome knives up to 185mm in length and
equipped with automatic timer to adjust sharpening time. It comes
with universal knife holder block, redressing pad, honing glass
plates, honing abrasive emery coarse and fine both and oil set

Cat No	Item	Type	Thin	Rs
AB-85-A	Automatic knife Razer sharpener for Microtome knife For Blade 185mm line with timer			
AB-85-S	Automatic knife Razer sharpener for Microtome knife Set honing glass plates, honing abrasive emery coarse and fine both and oil set			

