

A b

**Dissolved Oxygen DO Meter 0-20 ppm
With Probe**

With probe with membranes portable kit 9vDC
Complete set in portable case with manual

AC-303-A Dissolved Oxygen DO Meter Portable
with probe 9v battery and membranes in case



Abron Exports
www.abronexports.com
www.abronexports.com

Range	D.O	0 to 20 ppm
	Temp	0 to 50° C
Resolution	D.O	0.1 ppm
	Temp	0.1° C
Accuracy	D.O	0.1 ppm ± 1 count
Temp Compensation	0 to 50° C Manual	
Accuracy compensation	± 3%	
D.O Sensor	Amperometric Gold/Silver Membrane type	
Display LCD	3½ digit seven segment	
Power	1×9 V Battery	
Accessories	1) Gold/Silver D.O Probe 2) Dust Cover 3) Instruction Manual	

**Dissolved Oxygen DO Meter 0-20ppm x
LC0.1 Table Model With probe
With stand and membranes**

AC-303-B Dissolved Oxygen DO Meter 0.20ppm LC0.1
Table model with probe , membranes 220V



Range	0 to 20.0 ppm
Resolution	0.1 ppm
Accuracy	0.1 ppm ± 1 count
Temp compensation	0° to 50° C Manual
Accuracy of compensation D.O	± 3%
Display	3½ digit LED Display
Power	220V ± 10% Ac 50Hz
Accessories	1) Gold/Silver D.O Probe 2) Dust Cover 3) Instruction Manual

**Dissolved Oxygen chemical Test kit
With test tube etc with matching chart
AB-97DO Dissolved test kit chemicals set**

