

Karl Fisher Titrator Meter Digital Apparatus With Magnetic stirrer

Karl Fischer Titrator is an automatic apparatus designed to carry out Karl Fischer Titrations for accurate determination of moisture content in various samples having control and stirring with array of 10 LED's buzzer, auto zero burette, Adjustable timer & Dual Platinum electrode

Display	Array at 10 LED in line
Sample Size	1-50mg of water
K.F. Dispensing	Resolution 0.05ml
Indication	Run add End Display in Auto Mode
Detection	Current sensing by electronic control
K.F. Flow	Highly sensitive solenoid control
Electrode	Dual Platinum Electrode
Titration Cell	Air tight 200ml glass vessel
End Point Delay	30second
Buzzer	Yes
Timer	Adjustable Timer
Stirrer on/ Off	Variable speed magnetic stirrer
Burette	Pressure filling auto zero burette with 0.05ml resolution
Power	230V + 10% AC 50Hz

Accessories

- 1) Burette 10 ml
- 2) Reservoir Bottle 500 ml.
- 3) Rubber Bellow for Pumping Air
- 4) Titration Vessel 300 ml
- 5) Dual Platinum Electrode
- 6) Teflon Coated Magnetic Stirring Capsule
- 7) Moisture Trap with Rubber Tubing
- 9) Solenoid Valve
- 10) Instruction Manual

Cat No	Specifications	Type
AC-305-A	Digital Karl Fisher Titration + Auto Burettes With LED and controller 220v with bellow	



Kipps Gas Preparation Apparatus Boro Glass

Borosilicate glass flask special set with input funnel for acid input for solid iron sulphide, output H₂S gas by stop cock parts as per pic below

Cat No	Type	Borosilicate glass	Cap	Unit
AC-437-A	Gas Generator	Kipps apparatus	250 ml	
AC-437-B	Gas Generator	Kipps apparatus	500 ml	
AC-437-C	Gas Generator	Kipps apparatus	1000ml	
AC-437-D	Gas Generator	Kipps apparatus	2000ml	



Kipps Gas Preparation Apparatus Poly PP

polythene three main portions top fitted with a thistle funnel through which dilute acid reach Iron sulphide to base. H₂S gas collects in the middle it can be released through a leak proof stopcock self lubricant PTFE base has residual outlet for easy cleaning

Material : High Density Polyethylene

Cat No	Type	Borosilicate glass	Cap	Unit
AC-504-01	Gas Generator	Kipps apparatus	250 ml	
AC-504-02	Gas Generator	Kipps apparatus	500 ml	
AC-504-03	Gas Generator	Kipps apparatus	1000 ml	



Abron Exports
www.abronexports.com