

Microscope Inclined Monocular Research Pathological Bright field AB-10

High quality Inclined Head	Inclined Prismatic Tube
Stage 145x125mm	Rotatable Prismatic Tube through 360°
Condenser	Co-axial mechanical stage built in
Abbe Condenser	Sub stage movable by rack and pinion
Light control	1.25 N.A
Focusing Lockable	Iris Diaphragm and filter
Fine focusing	Coarse focusing by rack and pinion
Nose Piece	Graduated drum 0.002 mm.
Illumination	Quadruple revolving nose piece
Optics Eye Piece	light in base 6 volts, 20W variable
Objectives	10x WF & 15x Huygenian
	Achromatic 5x, 10x, 40 or 45x (SL)
	100x Oil Immersion spring loaded SL

Huygenian Eye Piece	Achromatic Objective	Magnification
10X WF	5X	50x 75X
Wide Field	10 X	100X 150 X
15X	45 X	450X , 675X
	100x Oil Immersion	1000x, 1500X

Cat No Type Range Optics

AB-10-A Microscope Monocular **Coaxial mechanical stage**
Iris diaphragm condenser lens focusable LED light
Eye piece 10x,15x Objective 10x,45x 100x Oil lens

With special Din Optics

AB-06-BD Medical Monocular Coaxial mechanical stage
Iris diaphragm condenser lens focusable LED Light
Din Optics Eye piece 10x,15x Objective 10x,45x 100x Oil lens



Microscope Inclined Monocular Research Pathological Bright field Microscope Special Coaxial body AB-10S

High quality Inclined Head	Inclined Prismatic Tube
Stage 145x125mm	Rotatable Prismatic Tube through 360°
Condenser	Co-axial mechanical stage built in
Abbe Condenser	Sub stage movable by rack and pinion
Light control	1.25 N.A
Focusing Lockable	Iris Diaphragm and filter
Fine focusing	Coarse focusing by rack and pinion
Nose Piece	Graduated drum 0.002 mm.
Illumination	Quadruple revolving nose piece
Optics Eye Piece	light in base 6 volts, 20W variable
Objectives	10x WF & 15x Huygenian
	Achromatic 5x, 10x, 40 or 45x (SL)
	100x Oil Immersion spring loaded SL

Huygenian Eye Piece	Achromatic Objective	Magnification
10X WF	5X	50x 75X
Wide Field	10 X	100X 150 X
15X	45 X	450X , 675X
	100x Oil Immersion	1000x, 1500X

Cat No Type Range Optics

AB-10-SA Microscope Monocular **Coaxial Body + stage**
Iris diaphragm condenser lens focusable LED light
Eye piece 10x,15x Objective 10x,45x 100x Oil lens

AB-10-SC Microscope Monocular **Coaxial Body + stage**
Iris diaphragm condenser lens focusable LED light
Eye piece 10x,15x Objective 10x,45x 100x Oil lens

With special Din Optics

AB-06-SBD Microscope Monocular **Coaxial Body + stage**
Iris diaphragm condenser lens focusable LED Light
Din Optics Semi plan optics
Eye piece 10x,15x Objective 10x,45x 100x Oil lens



Old AB-10SC



New AB-10SA

Microscope with Digital Camera Eye

Cat No with software for laptop Mega Pixel

AB-10D1 Microscope Monocular + Digital camera software 1.3MP
Eye piece 10x,15x Objective 10x,45x,100x Oil lens

AB-10D3 Microscope Monocular + **Digital camera software** 3 MP
Eye piece 10x,15x Objective 10x,45x,100x Oil lens

AB-10D5 Microscope Monocular +**Digital camera software** 5 MP
Eye piece 10x,15x Objective 10x,45x,100x Oil lens

Microscope with Digital Camera Eye

Cat No with software for laptop Mega Pixel

AB-10SD1 Microscope Monocular + Digital camera software 1.3MP
Eye piece 10x,15x Objective 10x,45x,100x Oil lens

AB-10SD3 Microscope Monocular + Digital camera software 3 MP
Eye piece 10x,15x Objective 10x,45x,100x Oil lens

AB-10SD5 Microscope Monocular +Digital camera software 5 MP
Eye piece 10x,15x Objective 10x,45x,100x Oil lens