

Millipore Filtration Sterility set for 47mm paper + Vacuum Pump setup Abron AC-471MS



Vacuum Filtration Manifold.

Materials:

- **Body and branches** : Type 316L stainless steel
- **Valves** : PTFE stopcock, Electroplated SS body.

Connections:

- **Filter Head** : 6 Stage, 3 Stage, Single Stage.
- **Supported Funnel** : Glass , PP And Stainless Steel.

Cat No	Type	Qty Set
AC-471-MS3	Filter holder manifold Multi	set of 3
AC-471-MS6	Filter holder manifold Multi	set of 6

<p>Bottle Vacuum wide mouth for Polypropylene with silicon cork and tube attachment pipes</p>  <p>Abron Exports www.abronexports.com</p>	<p>Borosilicate glass Conical flasks with silicon cork and tube</p>  <p>Abron Exports www.abronexports.com</p>	<p>Vacuum Tube Heavy 8 mm dia 2 meter Long</p>  <p>Abron Exports www.abronexports.com</p>
<p>Micro Filter paper filtration Membrane Cellulose Nitrate suitable for Millipore</p> <p>AC-607-B Filter Paper 2 μm Filtration Dia 47mm Pack of 100 AC-607-D Filter Paper 45μm Filtration Dia 47mm Pack of 100</p>  	<p>AC-288F Vacuum Pump Oil Free Rotary Diaphragm capacity 15 or 25 litre / Minutes</p> <p>Maximum Vacuum 554 mm of Hg - 22" Inch Maximum Pressure 25 Psig Motor 220v 50Hz/20 HP Approx weight 3 kg size mm 175x110x150mm</p> <p>AC-288F15 Cap 15 Liters / Minutes</p>  <p>AC-288F25 Cap 25 Liters / Minutes</p>	