

Micro Pipette Stand round 6 pipette Abron

With base and top with capacity 6 pipette

Cat No	Type	Capacity
AC-289-S6	Micro Pipette stand Round	6 pipette
AC-289-ST6	Micro Pipette stand Round SP	6 pipette



Micro Pipette Stand Abron (Z) Type

Perspex Z shape with notch 6 pipette

Cat No	Type	Capacity
AC-522-52	Micro Pipette Z perspex	6 Pipette
AC-289-SZ	Micro Pipette Z perspex	6 Pipette



Micro Pipette Stand Abron Round

Plastic L shape two colours with 6 pipette

Cat No	Type	Capacity
AC-289-SL	Micro Pipette L plastic	6 Pipette



Micro Pipette Tips General Abron

Compatible to all types of micro pipettes fixed/ variable

Cat No	Type	Size	colour	Pack
AC-532-PT-1	Pipette tips	0.2-10 ul	White	1000 Pc
AC-532-PT-2	Pipette tips	2-200 ul	Yellow	1000 Pc
AC-532-PT-3	Pipette tips	200-1000 ul	Blue	1000 Pc
AC-532-PT-4	Pipette tips	1-5 ml	White	100 Pc



Micro Tip Box For Pipette Tip

These boxes are molded in polypropylene and are therefore Autoclavable. The Box itself is clear providing an excellent view of the tips inside while a table like structure is provided inside to hold the tips. The tip holding spaces are uniformly placed to facilitate the use of multichannel pipettes.

Cat No	Item	Size	Cap box
AC-532-MC10	Pipette tip box	10 ul tip	96 Tips
AC-532-MC20	Pipette tip box	200 ul tip	96 Tips



Micro Tip Box For Pipette Tip

These boxes are molded in polypropylene and are therefore Autoclavable. The Box itself is clear providing an excellent view of the tips inside while a table like structure is provided inside to hold the tips. The tip holding spaces are uniformly placed to facilitate the use of multichannel pipettes.

Cat No	Item	Size	Cap box
AC-532-MC1H	Pipette tip box	1000 ul tip	36 Tips
AC-532-MCP	Pipette tip	Loose	100 Tips

