

# STEREO DISSECTING MICROSCOPE

## ABRON MODEL AB-33 TURRENT TYPE

### Fixed Objectives



### AB-33 RESEARCH STEREO BINOCULAR MICROSCOPE

- \* Turrent mount Achromatic Objective
- \* Precise accurate fine motion
- \* Top & Bottom built in lights with TXR
- \* One Eye Piece diopter adjustment with rubber guard
- \* Round glass stage changable with white perspex stage plate
- \* **DUEL LIGHT SYSTEM**  
Incident & trasmitted light  
6V bulb 220v operated.
- \* Thermocol packing With dust cover

#### AB-33B AS ABOVE BUT

Eyepieces WF 10x paired  
Objective 2X & 4X  
Magnification 20X & 40x

#### AB-33B AS ABOVE BUT

Eyepieces WF 10x paired  
Objective 3X & 5X  
Magnification 30X & 50x