


Millipore Filter Folder multi manifold SS For 47mm Micro filter paper

Abron AC-471MS3

Stainless steel holder with Teflon stopcock



Standard Configuration of Vacuum Filtration Manifold System



Cat No	Type	Qty Set
AC-471-MS3	Filter holder manifold Multi	set of 3
AC-471-MS6	Filter holder manifold Multi	set of 6

Bottle Vacuum wide mouth for Polypropylene with silicon cork and tube attachment pipes



Borosilicate glass Conical flasks with silicon cork and tube



Vacuum Tube Heavy 8 mm dia 2 meter Long



Micro Filter paper filtration Membrane Cellulose Nitrate suitable for Millipore

AC-607-B	Filter Paper 2 μ m Filtration Dia 47mm Pack of 100
AC-607-D	Filter Paper 45 μ m Filtration Dia 47mm Pack of 100



AC-288F Vacuum Pump Oil Free Rotary Diaphragm capacity 15 or 25 litre / Minutes

Maximum Vacuum 554 mm of Hg - 22" Inch
Maximum Pressure 25 Psig Motor 220v 50Hz/20 HP
Approx weight 3 kg size mm 175x110x150mm

AC-288F15 Cap 15 Liters / Minutes



AC-288F25 Cap 25 Liters / Minutes