

## Optical Bench Double Rod

### AP-840-A Optical Bench Double Bar 100cm or 150cm with 4 riders 2 pins steel graduated rod set

Two metallic rods (one engraved in cm mm) cast iron feet, accessories Parts are 4 metal riders Two fixed two having transverse motion by knob, two lens holders two optical pins

AP-840-B Optical Bench double rod 1 meter +parts


AP-840-BS Optical Bench double rod 1 meter +parts  
Both rods stainless steel



### AP-841-A Optical Bench Double Rod Heavy 100cm or 150cm with 4 riders 2 pins steel graduated rod set

1½ meter 150cm long double rod both rods covered with brass pipe np one with graduated in cm mm accessories parts heavy cast iron powder coated are 4 metal riders. Two fixed and two having transverse motion by knob, two lens holders, two optical pins thick all parts heavy accurate precise work  
AP-841BR Optical Bench double rod 150cm long  
With both rod covered with brass (1 graduated) pipe

AP-841-SS Optical Bench double rod 150cm long parts  
Both Rods (1 graduated ) of Stainless Steel metal



### AP-840- Optical Bench Other Options

AP-840-C Optical Bench double rod 1 meter + parts  
With one rod covered with Brass graduated pipe

AP-840-D Optical Bench double rod 150cm long

AP-840-E Optical Bench double rod 150cm long


With one rod covered with Brass graduated pipe

AP-840-F Optical Bench double rod 150cm long

With both rod covered with brass pipe

AP-840-G Optical Bench double rod 150cm long

Both Rods of Stainless Steel metal



### AP-841-H Optical Bench Extra Heavy with 3<sup>rd</sup> rod 150cm with 4 riders 2 pins steel graduated rod steel feet powder coated set

Two fixed two having transverse motion by knob

