



## Labomed Stereo Zoom Microscope CSM4 (Abron Cat AB-38B-4csm )

Modular stereo zoom Binocular and Trinocular microscope with a distinct combination of optical performance, flexibility and exceptional value. Designed with the Greenough optical system, CZM4 provides distortion free high resolution erect 3D images. The lateral zoom control adds convenience to operation. The choice in viewing-heads, stands, illumination and accessories accommodate a wide variety of operating conditions. At an economical price, the new CZM4 offers a basic platform for Zoom Stereoscopic format for industrial as well as educational needs. A host of accessories such as Auxiliary lenses, Digital camera attachment and a variety of stands make CZM4 the best choice of OEM requirements.

<b>Specifications</b>	
<b>STAND</b>	Dual light stable stand with focusing knobs. Optional stands with fluorescent light, Boom stand, with bonder arm, Pillar type stand, Slim Base Stand.
<b>VIEWING BODIES</b>	Binocular Zoom Pod 1:4, mag. 10x-40x, inclined at 45, IPD 54-74mm Trinocular Zoom Pod 1:4, mag. 10x-40x, inclined at 45, IPD 54-74mm
<b>ILLUMINATION</b>	Provided with built in illumination both for incident light(12V 20W Halogen) and sub-stage transmitted light(12V 10W).Optional Fluorescent Light.
<b>ZOOM RANGE</b>	1:4
<b>EYEPIECES</b>	Focusable WF 10x(FN 22mm). Optional WF 10x (FN22), WF 16x (FN16mm), WF 20x FN12mm) with fixed or focusable movement.
<b>AUXILIARY OBJECTIVES</b>	0.5X, 1.5X and 2.0X (Optional)
<b>ACCESSORIES</b>	Darkfield attachment, Polarising Kit, Ring Light, Fibre Optic Illumination, Dual light guides, photo adapters for 35mmSLR camera, Digital Still camera.

[Send Enquiry](#)