



AB-11C BINOCULAR RESEARCH MICROSCOPE CO-AXIAL STAGE LED

Focussing- for coarse/fine

Mechanical Stage-145x 125mm coaxial

Binocular Head-High Resoulution

Inclined Rotatable 360°

A built-in base LED with variable intensity

FILTER Blue and Green Filter

EYE PIECES- Wide field 10x paired

15x Hygenion. **PAIRED**

OBJECTIVE- 4x, 10x 40x,
100x Oil immersion

MAGNIFICATION - 40x, 60X

100x 150X 400x, 600X, 1000X, 1500X

