

Vacuum Desiccator 300mm + Vacuum Pump Oil Free 15 Lt Per Minute Gauge 760mm With Tube Set Abron AC-288K



AC-288M4B

+ AC-515PV

Abron Exports
www.abronexports.com
AC-288F15

Desiccator (Vacuum) with stop cock vacuum and perforated disc Polyethylene Polycarbonate AC-515PV

Autoclavable vacuum up to 740 mm of Hg 24 hrs
The top dome, moulded in rigid and transparent polycarbonate, gives a crystal clear view of the desiccant placed knob for easy handling lower part is of polypropylene. Silicon rubber O-ring Polypropylene Stopcock with leak proof PTFE plug Three way vacuum creations, shutting off & releasing

Material: Polyethylene/ Polycarbonate

Vacuum Pump Oil Free 15 Lit/ Min AC-288F15

Salient Features :

- Diaphragm Oil - Free Air Compressors
- Body Special Grade Aluminium
- Made Of 2 - Ply Nylon Reinforced Neoprene Rubber.
- Valves Are Ss 316
- Oil Free No Lubrication Is Required
- Noiseless Performance
- Diaphragms Teflon
- Maximum Air Flow 15 Litres/ Minutes
- Maximum Vacuum 554 Mm Of Hg - 22" Inch
- Maximum Pressure 25Psig Motor 220v 50hz 1/20HP
- Approx Weight 3 Kg Size Mm 175x110x150mm

<u>Cat No</u>	<u>Item</u>	<u>Type</u>	<u>Size</u>	<u>Oil Free Air Cap</u>
AC-288K30	Vacuum Desiccator	Dia 300mm	12 Inch	With Vacuum Pump Oil Free 15Lt/Min
AC-288K24	Vacuum Desiccator	Dia 240mm	10 Inch	With Vacuum Pump Oil Free 15Lt/Min
AC-288K20	Vacuum Desiccator	Dia 200mm	10 Inch	With Vacuum Pump Oil Free 15Lt/Min

OPTIONAL

AC-288M4B Vacuum Pump Oil Free Superior With Water Trap Gauge