

**Salinity Meter Digital LCD 1.999,19.99x 0.1
With Probe portable 9V battery AC-313SB**

Range 3 ranges 0 to 1.999 x 0.01 ppt
0 to 19.99 x 0.01 ppt
Display 3½ digit LCD display
Temp Compensation 0° C to 50° C manual
Cell Platinum Dip
Power 9 V battery
Accessories Salinity Cell with manual

Cat No	Item	Size	Thickness	Pack
AC-313-SB	Salinity meter Digital Portable 9V with probe	Range 1.999ppt, 19.99t x0.1ppt	abron	



**Salinity Meter Digital LED 1.999,19.99,49.99
LC 1digit with Probe 220v Table AC-313SC**

Range 3 ranges 0 to 1.999 x 0.01 ppt
0 to 19.99 x 0.01 ppt
0-49.99 x 0.01 ppt
Display 3½ digit LED display
Accuracy ± 2 FS ± 1 digit
Resolution 0.01 ppt
Temp Compensation 0° C to 50° C manual
Cell Constant Adjustable on Display
Cell Platinum Dip
Power 230 V ± 10% 50Hz
Dimensions 275x175x75 mm
Weight 2 Kg (approx.)
Accessories Salinity Cell and manual

Cat No	Item	Size	Thickness	Pack
AC-313-C	Salinity meter Digital table model manual with probe	Range 1.999ppt, 19.99t x0.1ppt , 0-50ppt x 1ppt	Manual temperature	



**Digital Conductivity TDS Salinity
Temperature ATC Meter with probes
Microprocessor Based Rs232 AP-301CP**

With ATC Auto Temperature probe

Conductivity

Range 0.1ms/cm 0 to 1000m mho in 6 ranges
Range selection automatic 6 ranges 0 to 1000ms/cm
Cell constant adjustable
Temperature co-off adjustable between 0 to 4%
Temp compensation automatic for entire range

TDS

Range 0 to 1000ppt in 6 ranges
Range selection automatic Resolution 0.1ppm

Salinity

Range 0 to 50ppt x 0.1 ppt
Accuracy ± 0.5% of range ± 1digit
Temp compensation Automatic
Salinity factor .75 and adjustable manually

Temperature

Range -20.0 to +200°c x 0.1°c
Input type PT 100 Calibration Manual
Display 16 character 2 lines alphanumeric LCD
Display format temp & conductivity or temp & TDS
temp & salinity simultaneously

Setting by soft touch keys
Storage 1000 samples with date & time
Interface RS-232, Centronics & computer
Error indication t-error, cal error
Print format a) All data b) data wise

Accessories

- 1) Conductivity Cell
- 2) Stand
- 3) Temp Probe
- 4) RS 232
- 5) Dust Cover
- 6) Instruction Manual

Cat No	Specifications	Range	Accuracy
AC-301-CP	Digital Salinity Temp Conductivity TDS LCD Meter with ATC automatic temperature TDS Probe With Printer Output RS232	16x2 line	

