

Senior Rotary Microtome 0-50 LC 1 Micron

Block Holder Blade L Block Abron AB-90SR



Cat No	Item	LC
AB-90SR	Microtome Senior Closed Pattern With Wheel To cut Sections Feeding 2-50 Micron LC 1 Micron 1 No Disposable Blade Holder 5 No Disposable Blade For Soft Tissue 3 No Object Holder 1 Set L Block Sides 2 No Block Making Mould SS 5 No Block Making Rings Plastic	

Microtome Optional spares AB-91K


Cat No	Item	Pack
AB-91KNL	Disposable blade low profile	50 Pc
AB-91KNH	Disposable blade high Profile	50 Pc



AB-91KNH Disposable blade Holder for Low/ High Profile



AB-91LM25
Microtome Block L moulds Brass



AB-91M Microtome SS wax Mould 100pc



AB-91MP Microtome PP Rings wax Mould 100pc

