

Hot Air Oven Thermostatic Controller
STAINLESS STEEL SHEET SS

Thermostat By Capillary Thermostat
Outer Body Mild steel Powder coated
Inner Chamber **STAINLESS STEEL SHEET**

Heaters NICHROME WIRE
Shelves Two Stainless steel Sheet
Operates 220V AC 50Hz

Temperature +5' RTM 50⁰C To 220⁰C

Accuracy ± 2% -1⁰C
Heat Switch 3 Way Rotary.
Air went Closable.
Cable with wire 2 meter and plug

Hot Air Oven Temperature Digital
Pid controlled STAINLESS STEEL SHEET SS
Digital Temperature Controller Cum Indicator

Outer Body Mild steel Powder coated
Inner Chamber **STAINLESS STEEL SHEET**

Heaters NICHROME WIRE
Shelves Two Stainless steel Sheet
Operates 220V AC 50Hz

Temperature +5' RTM 50⁰C To 220⁰C

Accuracy ± 2% -1⁰C
Heat Switch 3 Way Rotary.
Air went Closable.
Cable with wire 2 meter and plug

Cat No	Size (inch)	CM	Volume	Watts
AC-344S12	12"x12"x12"	300x300x300mm	28 litre	800W
AC-344S14	14"x14"x14"	350x350x350mm	45 litre	1000W
AC-344S16	16"x16"x16"	400x400x400mm	65 litre	1200W
AC-344S18	18"x18"x18"	450x450x450mm	95 litre	1500W
AC-344S21	21"x21"x21"	530x530x530mm	140 litre	1800W
AC-344S24	24"x24"x24"	600x600x600mm	224 litre	2000W
AC-344S24B	24"x24"x36"	600x600x900mm	336 litre	2500W

Cat No	Size (inch)	CM	Volume	Watts
AC-344SD12	12"x12"x12"	300x300x300mm	28 litre	800W
AC-344SD14	14"x14"x14"	350x350x350mm	45 litre	1000W
AC-344SD16	16"x16"x16"	400x400x400mm	65 litre	1200W
AC-344SD18	18"x18"x18"	450x450x450mm	95 litre	1500W
AC-344SD21	21"x21"x21"	530x530x530mm	140 litre	1800W
AC-344SD24	24"x24"x24"	600x600x600mm	224 litre	2000W
AC-344SD24	24"x24"x36"	600x600x900mm	336 litre	2000W
AC-344SD36	36"x36"x36"	900x900x1000mm	720 litre	2500W
AC-344SD40	40"x40"x100"	1000x1000x2500mm	2500 litre	3000W



OPTIONS

- AC-344-F Oven Thermostatic spare Air Fan
- AC-344-EH Oven Thermostatic spare Heavy Exhaust
- AC-344-T **Timer Dial 24H x15mt on off**
- AC-344-PID Digital PID Controller cum Indicator
- AC-344-PP Oven Special PID Controller Digital cum Indicator +Programmable
- AC-344-PD **Digital PID Controller cum Indicator + Programmable with data logger**

NOTE- Above can be supplied with capillary thermostat also

OPTIONS

- AC-344-F Oven Thermostatic spare Air Fan
- AC-344-EH Oven Thermostatic spare Heavy Exhaust
- AC-344-T **Timer Dial 24H x15mt on off**
- AC-344-PID Digital PID Controller cum Indicator
- AC-344-PP Oven Special PID Controller Digital cum Indicator +Programmable
- AC-344-PD **Digital PID Controller cum Indicator + Programmable with data logger**

NOTE- Above can be supplied with capillary thermostat also