

Stereo Zoom Binocular Microscope AB-34b2
Eye Piece 10X WF Ratio 1:6.4 Dual light
Reflected & Transmitted LED Illumination

Suitable For Precision Dissection / Inspection Specimens.

Standard Specifications

Anti fungal treated multiplayer coating optical used,
45 Degree inclined head 50 – 75mm comfort viewing,
One eyepiece tube dioptre adjustment for spectacle wearer

Binocular head	360° Turn able
Pole type simple stand	Angle of viewing:45°
Zoom ratio:	1 : 6.4
Eye Piece SWF	WF10X / 20mm
Objective:	0.7x-4.5x
Magnification	7x-45x
Inter pupillary Distance	Adjustment range:55-75mm
Working distance	100mm
Reflected transmitted	LED illumination Super wide field
Rack & Pinion focusing	large Dia knobs
Stage	95mm Dia glass stage
Black & white	Perspex plate for field contrast

Accessories, Dust cover, manual, Thermocole box

AB-34-A2 Delux

Eye Piece WF	Objective Zoom	Magnification
10X Paired	1 : 6.4	
	0.7X to 4.5X	7X to 45X

Cat No	Type	Zoom	Ration
AB-34-A2	Stereo Zoom Binocular Microscope Ratio	1:6.4	
	Eye Piece 10x Pair Magnification Obj	0.7x-4.5X	
	With Dual Variable light in base LED/ 6V bulb		



Cat No	Type	Zoom	Ration
AB-34-T	Stereo Zoom Trinocular Microscope Ratio	1:6.4	
	Eye Piece 10x Pair Magnification Obj	0.7x-4.5X	
	With Dual Variable light in base 6V bulb / LED		

Stereo Zoom Trinocular Digital Microscope
With Digital Camera Eye AB-34D

Digital Eye Camera (remove eye piece and put it)
Or use in Trino Tube with camera + CD Software USB Cable
Can be used on Windows and Mac Book Apple in VLC Player

Binocular head	360° Turn able
Pole type simple stand	Angle of viewing:45°
Zoom ratio:	1 : 6.4
Eye Piece SWF	WF10X / 20mm
Objective:	0.7x-4.5x
Magnification	7x-45x
Inter pupillary Distance	Adjustment range:55-75mm
Working distance	100mm
Reflected transmitted	LED illumination Super wide field
Rack & Pinion focusing	large Dia knobs
Stage	95mm Dia glass stage
Black & white	Perspex plate for field contrast

Accessories, Dust cover, manual, Thermocole box

Dual position light	sliding system
Light 70% light to camera	30% to Binocular eyepiece

AB-34-D Digital Trino Models

Eye Piece WF	Objective Zoom	Magnification
10X Paired	1 : 6.4	
	0.7X to 4.5X	7X to 45X

Cat No	Type	Zoom	Ration
AB-34-D	Stereo Zoom Trinocular Microscope Ratio	1:6.4	
	Eye Piece 10x Pair Magnification Obj	0.7x-4.5X	
	With Dual Variable light in base LED/ 6VBulb		
AB-34-D3	Stereo Zoom Trini Microscope +Digital Camera		
	AB-14a 3 mp + software CD for laptop		
AB-34-D5	Stereo Zoom Trini Microscope +Digital Camera		
	AB-14a 5 mp + software CD for laptop		
AB-34-D10	Stereo Zoom Trini Microscope +Digital Camera		
	AB-14a 10 mp + software CD for laptop		
AB-34-D14	Stereo Zoom Trini Microscope +Digital Camera		
	AB-14a 14 mp + software CD for laptop		

Snap Shot Digital Camera Options

AB-34-DN	Trino Micro +camera digital snapshot	20MP NIKON
AB-13-DS	Trino Micro + camera digital snapshot	20MP SONY

Optional Digital DSLR with special tub can be supplied



Note :- laptop is not standard part of supply (optional)