



FAD4122200	ZOOM STEREO MICROSCOPE (BIOLOGY & INDUSTRIAL) CZM4	
	<p>CZM4 BINOCULAR ZOOM STEREO MICROSCOPE</p> <p>Stand: Dual light stable stand with Ball bearing focusing controls</p> <p>Viewing Bodies: Binocular tube, 45° inclined, 360° rotatable, IPD 54-74mm</p> <p>Eyepieces: Wide Field eyepiece WF 10x/22mm focusable paired with eye guards</p> <p>Magnification: 10x-40x</p> <p>Illumination: Built-in illumination system 12V-20W Halogen incident light and transmitted 12V-10W halogen transmitted light.</p> <p>ElectricalPower supply (220V-240V). Packed in styrofoam box, with operation manual, allen wrench, dust cover, cleaning cloth, power cord and spare fuses</p> <p></p>	
FAD4122200T	Same as above but with Trinocular observation head	

Magnification Table							
Objectives	Working Distance (in mm)	Eyepieces					
		W.F. 10x/18		W.F. 10x/22		W.F. 16x/16	
		Ob. Field (mm)	Mag.	Ob. Field (mm)	Mag.	Ob. Field (mm)	Mag.
None	90.0	5.5-22	40-10	4-16	64-16	3-12	80-20
0.5x	135.0	11-44	20-5	8-32	32-8	6-24	40-10
1.5x	43.0	3.7-14.7	60-15	2.7-10.7	96-24	2-8	120-30
2x	35.0	2.75-11	80-20	2.0-8.0	128-32	1.5-6	160-40

ABRON EXPORTS

E-mail No:- info@abronexports.com website <http://abronexports.com>