

A b

**Dissolved Oxygen DO Meter 0-20ppm x LC0.1
ATC Auto Temperature Model DO probe
With stand and Membranes AC-303C**

AC-303-C Dissolved Oxygen DO Meter 0.20ppm LC0.1
ATC Table model with probe , membranes 220V



Range	0 to 20.0 ppm
Resolution	0.1 ppm
Accuracy	0.1 ppm \pm 1 count
Temp compensation	0° to 50° C Manual ATC Probe
Accuracy of compensation D.O	\pm 3%
Temperature Range	0° to 50° C
Resolution	0.1° C
Accuracy	0.20% \pm 1 Count
Temperature Sensor	RTD Probe
Display	3½ digit LED Display
Power	220V \pm 10% Ac 50Hz
Accessories	1) Gold/Silver D.O Probe 2) Temperature Probe 3) Dust Cover 4) Instruction Manual .


**AC-303P
AC-303M**

**DO meter spare probe with 5 membrane
DO meter spare membranes**



**Dissolved Oxygen DO Meter 0-20ppm x 0.1
ATC -20-100c With DO probe+ stand
Membranes Printer RS232 computer Out
Microprocessor based RS232 AC-303M**

AC-303-M Dissolved Oxygen DO Meter 0.20ppm LC0.1
ATC Table model with probe membranes 220V
Microprocessor based



Range	0 to 20.0 ppm LC 0.01ppm
Temp compensation	Automatic & Manual With Keys
Pressure Adj	0 to 760 mm of Hg Automatic
Salinity Correction	0 to 20gm/l (Manual)
Calibration	Automatic
Temperature	Input PT100
Range	-20° to 100° Cx0.1°C
Display	16 Character * 2 Lines Alphanumeric LCD
Display Format	Temperature And DO both
Setting	By Soft Touch Keys
Storage	1000 Samples With Date & Time
Interface	RS-232, Centronics Computer
Real Time Display	24 Hour Mode With Date
Print format	a) All data b) Date Wise
Power Supply	220V \pm 10% Ac 50Hz
Accessories	1) Gold/Silver D.O Probe 2) Temperature Prob 3) Dust Cover 4) Instruction Manual .