



**AB-12D CAMERA FOR +BINOCULAR RESEARCH MICROSCOPE**

- +45mm DIN Achromatic Objective:4X,10X,40X(S),100X(S, Oil)
- + Wide Field Plane-Scope Eyepiece:WF10X
- +Total Magnification:40X-1000X
- +Articulated Free Binocular Head,
- +Interpupillary Distance:55-75mm,30° Inclined,360° Rotatable
- +Double Layers Mechanical Stage:125X125mm, Moving Range X-Y:60X40mm
- +Ø2-Ø30mm Iris Diaphragm
- +Ø32 Filter,ABBE Condenser:N.A. 1.25
- +Coaxial Coarse And Fine Focusing Adjustable Mechanism:14mm
- +Built-in Adjustable Brightness High intensity white light LED



AB-12DA CAMERA 1.5MP + software CD for laptop

AB-12DB CAMERA 3MP + software CD for laptop

AB-12DC CAMERA 5MP + software CD for laptop

AB-12Dd CAMERA 5MP + software CD for laptop  
High resolution

