



| | | |
|-----------------|-------------|-------------------|
| PART NO. | CZM6 | List Price |
|-----------------|-------------|-------------------|

| | | |
|--------------------------|--|--|
| FAD4122200 B6 | ZOOM STEREO MICROSCOPE (BIOLOGY & INDUSTRIAL) | |
| | <p>Stand: Dual light stable stand with Ball bearing focusing controls</p> <p>Viewing Bodies: Binocular tube, 45° inclined, 360° rotatable, IPD 54-74mm</p> <p>Eyepieces: Wide Field eyepiece WF 10x/22mm focusable paired with eye guards</p> <p>Zoom ratio: 1:6.9</p> <p>Magnification: 6.5x-45x</p> <p>Illumination: Built-in illumination system 12V-20W Halogen incident light and transmitted 12V-10W halogen transmitted light.</p> <p>Electrical: Power supply (220V-240V).</p> <p>Packed in Styrofoam box, with operation manual, allen wrench, dust cover, cleaning cloth, power cord and spare fuses</p> <p style="text-align: center;"></p> |  <p style="text-align: center;">Price Rs. 44460.00</p> |
| FAD4122200 T6 | Same as above but with Trinocular observation head | 49590.00 |

| Magnification Table | | | | | | | |
|---------------------|--------------------------|----------------|----------|----------------|----------|----------------|----------|
| Objectives | Working Distance (in mm) | Eyepieces | | | | | |
| | | W.F. 10x/22 | | W.F. 16x/16 | | W.F. 20x/12 | |
| | | Ob. Field (mm) | Mag. | Ob. Field (mm) | Mag. | Ob. Field (mm) | Mag. |
| None | 100.0 | 4.4-30.7 | 45-6.5 | 3.5-24.6 | 72-10.4 | 2.9-2.0 | 90-13 |
| 0.5x | 160.0 | 8.9-61.4 | 22.5-3.3 | 7.1-49.2 | 36-5.2 | 5.8-40 | 45-6.5 |
| 1.5x | 40.0 | 3.0-20.5 | 67.5-9.8 | 2.4-16.4 | 108-15.6 | 1.9-13.3 | 135-19.5 |
| 2x | 25.0 | 2.2-15.4 | 90-13 | 1.8-12.3 | 144-20.8 | 1.4-10 | 180-26 |

Telephone No. 0091 171 2699478

Fax No: 0091 171 2698459

ABRON EXPORTS

New Add. :- 105 Industrial Estate HSIIDC, AMBALA CANTT. 133 001

E-mail No:- info@abronexports.com website <http://abronexports.com>