

**Dissecting Microscope Bakelite Parts Body Die cast Eye piece 10x 20x with side rest AB-23A**

|                  |                                   |
|------------------|-----------------------------------|
| Body             | Cast Iron round base              |
| Bakelite parts   | Yes                               |
| Stage 4" x 4"    | glass with removable white plate  |
| Focussing        | By rack and pinion                |
| Illumination     | mirrors (one plane & one concave) |
| Optics Eye piece | 10x, 20x.                         |

| Cat No  | Type                                      | Parts | Eye piece |
|---------|-------------------------------------------|-------|-----------|
| AB-23-A | Dissecting Microscope Rack & pinion Brass |       | 10x 20x   |



**Dissecting Microscope Brass Parts Die cast Eye piece 10x 20x with side rest AB-23b**

|                  |                                   |
|------------------|-----------------------------------|
| Body             | Cast Iron round base              |
| Brass parts      | Yes                               |
| Stage 4" x 4"    | glass with removable white plate  |
| Focussing        | By rack and pinion                |
| Illumination     | mirrors (one plane & one concave) |
| Optics Eye piece | 10x 20x                           |

| Cat No  | Type                                      | Parts | Eye piece |
|---------|-------------------------------------------|-------|-----------|
| AB-23-B | Dissecting Microscope Rack & pinion Brass |       | 10x 20x   |



**Dissecting Microscope Brass Parts ISI Die cast Eye piece 10x 20x with side rest AB-24**

|                  |             |                                   |
|------------------|-------------|-----------------------------------|
| Body             | Brass Parts | Cast Iron round base              |
| Stage 4" x 4"    |             | glass with removable white plate  |
| Focussing        |             | By rack and pinion                |
| Illumination     |             | mirrors (one plane & one concave) |
| Optics Eye piece |             | 10x 20x                           |

| Cat No  | Type                                          | Parts | Eye piece |
|---------|-----------------------------------------------|-------|-----------|
| AB-23-B | Dissecting Microscope Rack & pinion ISI Brass |       | 10x 20x   |



**Dissecting Microscope Bull Lens Brass Rack & Pinion Parts AB-25**

|                  |                                   |
|------------------|-----------------------------------|
| Body             | Cast Iron round base              |
| Brass parts      | Yes                               |
| Stage 4" x 4"    | glass with removable white plate  |
| Focussing        | By rack and pinion                |
| Illumination     | mirrors (one plane & one concave) |
| Optics Bull Lens | Eye piece 83mm 10X                |

| Cat No | Type                                               | Parts | Eye piece |
|--------|----------------------------------------------------|-------|-----------|
| AB-25  | Dissecting Microscope Rack & pinion Bull lens 83mm |       | 10x       |



**Dissecting Microscope High quality Brass Part Die cast Eye piece 10x 20x +side rest Abron AB-26**

|                     |                                   |
|---------------------|-----------------------------------|
| Body                | Cast Iron round base              |
| Brass parts         | Yes                               |
| Stage 4" x 4"       | glass with removable white plate  |
| Focussing           | By rack and pinion                |
| Illumination        | mirrors (one plane & one concave) |
| Optics High quality | Eye piece 10x and 20X             |

| Cat No | Type | Parts | Eye piece |
|--------|------|-------|-----------|
|--------|------|-------|-----------|

AB-27 Dissecting Microscope Rack & pinion High lens 10x 20x



Abron Exports  
www.abronexports.com



www.abronexports.com